




SORATH ORTHO

A DIVISION OF SIGMA SURGICAL

Twin Lock-BCP Medial Distal Tibial Plate

Twin Lock - Bone Compression Plate System




ISO 13485

www.sorathortho.co.in

DISTAL TIBIA

Twin Lock-BCP Medial Distal Tibial Plate



Primary use : Distal Tibia Fixation
Screw Size : Twin Lock-BCP Screw \varnothing 3.5mm/ Twin Lock-BCP Cancellous Screw \varnothing 3.5mm / Cortical Bone Screw \varnothing 3.5mm/ Twin Lock-BCP Screw, \varnothing 3.7 mm, Cannulated
Matching Instrument : J0200016 Complete Instrument Set for 3.5mm Locking System

Size (Holes)	Reference Number		Size (Holes)	Reference Number	
	Titanium	Stainless Steel		Titanium	Stainless Steel
	Left			Right	
4	C01114004	C01111004	4	C01124004	C01121004
6	C01114006	C01111006	6	C01124006	C01121006
8	C01114008	C01111008	8	C01124008	C01121008
10	C01114010	C01111010	10	C01124010	C01121010
12	C01114012	C01111012	12	C01124012	C01121012
14	C01114014	C01111014	14	C01124014	C01121014

end to end

GLOBAL MEDICAL SUPPLY SOLUTIONS

○ Sourcing ○ Production ○ Quality Control ○ Logistics ○ Customer Services



- * Hyper, Sorath, Lexa, Spyrox, Sigma, Pelviflex and Flexora are trademarks or registered trademarks of SIGMA SURGICAL PVT. LTD. in INDIA, the USA and other Countries.
- * The information in the Catalogue is valid as of OCTOBER 2024 Designs and specifications are subject to be changed without notice.
- * Product shown in Catalogue may differ from actual products.
- * Sigma Surgical is not responsible for differences between the information in the Catalogue and the actual implants.

Manufacturing Facility :

SIGMA SURGICAL PRIVATE LIMITED

Survey no:1636/1,Opp. Bharat petrol pump,
Indore highway, Village-Kuha,
Taluka-Daskroi, Ahmedabad-382433.
Gujarat, India. 📍 : 22.99946, 72.78633

Customer care :

+91-7575011006 / +91-8866211414

www.sorathortho.com
info@sorathortho.com