




SORATH ORTHO

A DIVISION OF SIGMA SURGICAL

Twin Lock-BCP Proximal Medial Tibia Plate, Ø 3.5mm

Twin Lock - Bone Compression Plate System



ISO 13485

www.sorathortho.co.in

PROXIMAL TIBIA

Twin Lock-BCP Proximal Medial Tibia Plate, Ø 3.5mm



- Primary use** : Proximal Tibia Fixation
Screw Size : Twin Lock-BCP Screw Ø3.5mm/ Twin Lock-BCP Cancellous Screw Ø3.5mm / Cortical Bone Screw Ø3.5mm/
 Twin Lock-BCP Screw, Ø 3.7 mm, Cannulated
Matching Instrument : J0200016 Complete Instrument Set for 3.5mm Locking System

Size (Holes)	Reference Number		Size (Holes)	Reference Number	
	Titanium	Stainless Steel		Titanium	Stainless Steel
	Left			Right	
4	C02814004	C02811004	4	C02824004	C02821004
6	C02814006	C02811006	6	C02824006	C02821006
8	C02814008	C02811008	8	C02824008	C02821008
10	C02814010	C02811010	10	C02824010	C02821010
12	C02814012	C02811012	12	C02824012	C02821012
14	C02814014	C02811014	14	C02824014	C02821014

end to end

GLOBAL MEDICAL SUPPLY SOLUTIONS

○ Sourcing ○ Production ○ Quality Control ○ Logistics ○ Customer Services



- * Hyper, Sorath, Lexa, Spyrox, Sigma, Pelviflex and Flexora are trademarks or registered trademarks of SIGMA SURGICAL PVT. LTD. in INDIA, the USA and other Countries.
- * The information in the Catalogue is valid as of OCTOBER 2024 Designs and specifications are subject to be changed without notice.
- * Product shown in Catalogue may differ from actual products.
- * Sigma Surgical is not responsible for differences between the information in the Catalogue and the actual implants.

Manufacturing Facility :

SIGMA SURGICAL PRIVATE LIMITED

Survey no:1636/1,Opp. Bharat petrol pump,
Indore highway, Village-Kuha,
Taluka-Daskroi, Ahmedabad-382433.
Gujarat, India. 📍 : 22.99946, 72.78633

Customer care :

+91-7575011006 / +91-8866211414

www.sorathortho.com
info@sorathortho.com