


## Twin Lock-BCP Anarel Plate, Ø3.5mm Twin Lock - Bone Compression Plate System





www.sorathortho.co.in

### ANKLE

### Twin Lock-BCP Anarel Plate, Ø3.5mm



Primary use Screw Size

: Calcaneus Fixation

: Twin Lock-BCP Screw Ø3.5mm / Twin Lock-BCP Cancellous Screw Ø3.5mm Matching Instrument : J0200016 Complete Instrument Set for 3.5mm Locking System

Size (Holes)	Reference Number		Size	Reference Number	
	Titanium	Stainless Steel	(Holes)	Titanium	Stainless Steel
	Left			Right	
64	C02015064	C02011064	64	C02025064	C02025064
69	C02015069	C02011069	69	C02025069	C02025069
76	C02015076	C02011076	76	C02025076	C02025076
81	C02015081	C02011081	81	C02025081	C02025081
	I		I	I	

# end to end **GLOBAL MEDICAL SUPPLY SOLUTIONS**

OSourcing OProduction OQuality Control OLogistics OCustomer Services



\* Hyper, Sorath, Lexa, Spyrox, Sigma, Pelviflex and Flexora are trademarks or registered trademarks of SIGMA SURGICAL PVT. LTD. in INDIA, the USA and other Countries.

- The information in the Catalogue is valid as of OCTOBER 2024 Designs and specifications are subject to be changed without notice
  Product shown in Catalogue may differ from actual products.
  Sigma Surgical is not responsible for differences between the information in the Catalogue and the actual implants.

Manufacturing Facility :

#### SIGE SURGICAL PRIVATE LIMITED

Survey no:1636/1,Opp. Bharat petrol pump, Indore highway, Village-Kuha, Taluka-Daskroi, Ahmedabad-382433. Gujarat, India. **9** : 22.99946, 72.78633

Customer care : +91-7575011006 / +91-8866211414

www.sorathortho.com info@sorathortho.com