


Twin Lock-BCP Xyston Distal Humerus Plate Twin Lock - Bone Compression Plate System





(€ ISO 13485

DISTAL HUMERUS

Twin Lock-BCP Xyston Distal Humerus Plate



Primary use : Distal Humerus Fixation

Screw Size : Twin Lock-BCP Screw Ø3.5mm/ Twin Lock-BCP Cancellous Screw Ø3.5mm / Cortical Bone Screw Ø3.5mm/

Twin Lock-BCP Screw, Ø 2.7 mm

Matching Instrument : J0200019 Complete Instrument Set for Twin Lock-bcp Distal Humeral Plates

Size (Holes)	Reference Number		Size	Reference Number	
	Titanium	Stainless Steel	(Holes)	Titanium	Stainless Steel
	Left			Right	
3	C01314003	C01311003	3	C01324003	C01321003
5	C01314005	C01311005	5	C01324005	C01321005
7	C01314007	C01311007	7	C01324007	C01321007
9	C01314009	C01311009	9	C01324009	C01321009
11	C01314011	C01311011	11	C01324011	C01321011
	Left with Lateral Support			Right with Lateral Support	
3	C01214003	C01211003	3	C01224003	C01221003
5	C01214005	C01211005	5	C01224005	C01221005
7	C01214007	C01211007	7	C01224007	C01221007
9	C01214009	C01211009	9	C01224009	C01221009
11	C01214011	C01211011	11	C01224011	C01221011

end to end

GLOBAL MEDICAL SUPPLY SOLUTIONS

OSourcing OProduction OQuality Control OLogistics OCustomer Services



- * Hyper, Sorath, Lexa, Spyrox, Sigma, Pelviflex and Flexora are trademarks or registered trademarks of SIGMA SURGICAL PVT. LTD. in INDIA, the USA and other Countries.
- * The information in the Catalogue is valid as of OCTOBER 2024 Designs and specifications are subject to be changed without notice
 Product shown in Catalogue may differ from actual products.

 Sigma Surgical is not responsible for differences between the information in the Catalogue and the actual implants.

Manufacturing Facility:

SIGE SURGICAL PRIVATE LIMITED

Survey no:1636/1,Opp. Bharat petrol pump, Indore highway, Village-Kuha, Taluka-Daskroi, Ahmedabad-382433. Gujarat, India. ♥: 22.99946, 72.78633

Customer care : +91-7575011006 / +91-8866211414

www.sorathortho.com info@sorathortho.com